

RT EPID-2

Radiotherapy Phantom

RT EPID-2 Tests

- Planar MV Imaging QA
- Geometric accuracy
- Low contrast detectability
- Uniformity and noise





RT EPID-2

Radiotherapy Phantom

RT EPID-2 Features

- 220mm x 220mm x 24mm
- Alignment marks for room lasers
- 6x 2mm Tungsten sphere, 2x 1mm Tungsten spheres for point spread function (PSF) and geometric accuracy test
- 4x ROIs of Aluminium for uniformity test
- 30 low contrast resolution details of diameter 1.0, 2.0, 4.0, 7.0, 10.0, 15.0 mm

